

CLAIMS

1. A quality assurance system for the detection of proliferative microorganisms, comprising
 - 5 a) a system for enriching microorganisms in a sample in an "overnight culture" corresponding to 8 to 24 hours' cultivation under standard conditions according to International pharmacopoeias (for example Ph. Eur.), food laws, cosmetics directives or other commercially available indirect methods,
 - 10 b) a kit for detecting living, damaged or dead microorganisms in filterable and/or non-filterable products containing
 - i) at least one reagent containing an inductor and a fluorescent reagent which, with living cells, leads to the formation of a certain enzyme that releases a detectable fluorescent dye by
15 reaction with a specific fluorescent reagent,
 - ii) at least one nucleic acid probe for detecting microorganisms by in situ hybridization, the nucleic acid probe being fixed to a fluorescent marker,
- 20 in which a detection limit for proliferative microorganisms of < 10 CFU/g is achieved.
2. A quality assurance system as claimed in claim 1, characterized in that a detection limit for proliferative microorganisms of < 1 CFU/g is achieved.
3. A kit according to claim 1, characterized in that the reagent from b i)
25 can be used for filterable liquid samples and products or for filterable liquid parts of the samples and products to be analyzed for detecting living microorganisms and indirectly for detecting dead microorganisms.
4. A kit according to claim 1, characterized in that the nucleic acid probe from b ii) can be used both for filterable liquid samples and products
30 and for non-filterable samples and products and for mixtures of filterable and non-filterable samples and products for detecting living

microorganisms.

5. A quality assurance system as claimed in any of claims 1 to 4 for detecting gram-positive and/or gram-negative bacteria and/or yeasts and/or molds and/or algae.

5 6. The use of the quality assurance system claimed in any of claims 1 to 5 for detecting microorganisms and for the quality assessment of filterable and/or non-filterable products and for evaluating the hygiene status of production plants, a detection limit of < 10 CFU/g being achieved.

7. The use of the quality assurance system claimed in any of claims 1
10 to 5 for detecting microorganisms and for the quality assessment of filterable and/or non-filterable products selected from the group consisting of crude products, cosmetic products, pharmaceutical preparations, foods, food supplements, beverages, textile auxiliaries, detergents and dyes and lacquers.

15 8. A process for detecting microorganisms in filterable and/or non-filterable products, in which the quality assurance system claimed in any of claims 1 to 5 is used by

a) cultivating the samples in an "overnight culture" corresponding to 8
20 to 24 hours' cultivation under standard conditions according to International pharmacopoeias (for example Ph. Eur.), food laws, cosmetics directives or commercially available indirect methods in order to enrich microorganisms and

b) using a kit for detecting living, damaged or dead microorganisms in
25 filterable and/or non-filterable products by

i) incubating the enriched sample with a reagent containing an
inductor and a fluorescent reagent which induces the
formation of a certain enzyme in the cells and, in the process,
allows a fluorescent compound to be formed from a
30 fluorescent reagent and/or

- ii) after fixing the bacteria, incubating them with a nucleic acid probe which is provided with a fluorescent marker in order to induce hybridization and
 - c) detecting the fluorescence of the samples and correlating the result with the number of cells, the number of cells being determinable and, where b) i) is used, dead and living cells being distinguishable from one another.
- 5